

# Original protocol used in Covassin et al. (2006) Dev. Biol. 299:551-562

## **Cell disruption for flow sorting modified from Zebrafish Book protocol**

### **You will need:**

Fetal Calf Serum  
1M CaCl<sub>2</sub>

### **Calcium free Ringer' solution**

116 mM NaCl  
2.9 mM KCl  
5 mM HEPES, pH 7.2.

### **Suspension medium**

Colorless (no phenol red) Leibovitz medium L-15, contains 0.3 mM glutamine  
(from GIBCO)  
0.8 mM CaCl<sub>2</sub>  
Pen 50 U/mL / Strep 0.05mg/mL (our stock is 10X)  
1 % FCS

### **Culture medium**

L15  
0.3 mM glutamine  
0.8 mM CaCl<sub>2</sub>  
Pen/Strep

### **Embryo Extract**

- Chill 200 3 day embryos after removing from chorions.
- Rinse in 0.5% chilled bleach for 2 min and then in zero calcium Ringer for 2 min.
- Transfer to a Dounce homogenizer with a minimum of liquid and homogenize well.
- Resuspend in 1 ml Culture Medium
- Store at -20°C

### **Protease solution**

0.25% trypsin (trypsin we have is 10X, 2.5%)  
1mM EDTA  
in PBS, pH=8

### **Collection medium**

Culture medium  
10% FCS  
10% Embryo Extract

### **Trizol (Invitrogen)**

**Do not forget the GFP- control (GFP- fish)!!**

**Protocol**

Grow embryos in egg water.

Dechorionate with pronase.

When embryos reach desired developmental stage, transfer them in 1.5mL tube and rinse them for 15min in calcium free Ringer.

**During this time**, get rid of yolk:

Cut very end of 200uL tip and make 150uL aliquots in 1.5mL tubes; for each do 3-5 gentle up and downs (make sure all the yolks are gone). Pool embryos, rinse well with Ringer.

Transfer into small culture dish with 5-10mL of Protease solution. Incubate at 28C. Every 10min homogenize well doing up and downs with 200uL tip. Monitor dissociation with microscope; it takes 30-50min.

Add 1-2mM CaCl<sub>2</sub> and 5-10%FCS to stop reaction

Centrifuge 3min at 3000rpm.

Discard supernatant and rinse cells in 10mL suspension solution. Centrifuge again.

Discard sup and resuspend cells in suspension medium  $10 \times 10^6$  cells/mL. At 24hpf, we approximate 25,000 cells/embryo.

For the control, if less than  $10 \times 10^6$  cells, resuspend in 1mL.

Bring cells to flow sorting facility (5<sup>th</sup> floor Medical School Building, out of the elevator right and right), keeping them warm.

Bring along, some suspension buffer, in case she needs to dilute cells, and collection buffer (5mL/sample).

The experiment is done sterile at room temperature. She collects all GFP+ and in some cases she has to collect GFP- as well (at least as much as GFP+).

When you get the samples back, centrifuge, discard sup, resuspend cells in 250uL of Trizol and transfer in 1.5mL tube. Store at -80C.

Proceed to total RNA isolation.